

Premium Restorations

	Premium PFM	Pressed e.max™	e.max CAD™	Crystalite	Zirconia (Layered)	Lava™ Esthetic
Material	Semi-precious or precious alloy coping with layered porcelain overlay	Ivoclar™ e.max lithium disilicate	Ivoclar™ e.max lithium disilicate	Monolithic super translucent zirconia	Zirconia (layered) coping with porcelain overlay	Esthetic Fluorescent Full-Contour Monolithic Zirconia
Method of manufacture	Traditionally-cast & built up by senior technician	Waxed & pressed by hand	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed & manufactured (Framework) & built up by hand	Digitally designed and manufactured
Aesthetics	★★★★☆	★★★★★	★★★★★	★★★★★	★★★★☆	★★★★☆
Shades / colouring <small>(using Vita Classical, 3D-Master or Chromascop guides)</small>	Internal effects customised shading multi-layered	Customised shading	Customised shading	No internal effects standard shades	No internal effects standard shades	No internal effects, standard shades
Strength	★★★★☆☆	★★★★☆☆	★★★★☆☆	★★★★★	★★★★☆☆	★★★★☆☆
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	✓	✗	✗	Fit to denture	Fit to denture, rests and ledges	Fit to denture
Benefits	Traditional method choice of alloys 10x microscope fit multi-layer build	Premium aesthetics metal-free	Premium aesthetics excellent fit metal-free	High durability extreme strength metal-free, exceptional accuracy	Good fit metal-free	High durability strength metal-free, exceptional accuracy
Recommended use	Full mouth, crowns & bridges	Full mouth, crowns, inlays/onlays veneers, 3/4 crowns	Full mouth, crowns, inlays/onlays veneers, 3/4 crowns	Full mouth, single crowns	Full mouth, crowns & bridges	Full mouth, single crowns
Warranty	6 years	1 year	5 years	15 years	1 year	10 years

Standard Restorations

	<p>Standard PFM</p>	<p>Full Gold Crown</p>	<p>Opalite</p>	<p>Lava™ Plus</p>	<p>Enamic</p>
Material	Semi-precious or precious alloy coping with porcelain overlay	Non-precious (gold-coloured) or precious (gold) alloy CoCr (silver coloured)	Monolithic translucent zirconia	Monolithic translucent zirconia	Hybrid Ceramic
Method of manufacture	Traditionally-cast & built up	Traditionally-cast	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed and manufactured
Aesthetics	★★★★☆	★☆☆☆☆	★★★★★	★★★★★	★★★★☆
Shades / colouring <small>(using Vita Classical, 3D-Master or Chromascop guides)</small>	No internal effects standard shades	Not required	No internal effects standard shades	No internal effects standard shades	No internal effects standard shades
Strength	★★★★☆	★★★★★★	★★★★★★	★★★★★★	★★★★★
Add-on options <small>(fit to denture, attachments etc. Charges apply)</small>	✓	✓	Fit to denture	Fit to denture	✗
Benefits	Traditional method choice of alloys 5x microscope fit	5x microscope fit exceptional strength	High durability extreme strength metal-free, exceptional accuracy	High durability extreme strength metal-free, exceptional accuracy	Minimal invasive. shock absorbing
Recommended use	Full mouth, crowns & bridges	Posterior teeth	Full mouth, crowns & bridges	Full mouth, crowns & bridges	Full mouth, crowns, inlays, onlays, implants
Warranty	6 years	6 years	15 years	1 year	5 years