







Splint Comparison Chart

						
	Occlusal Splints (Heat Cured)	Hard Soft	Soft	Anterior Mini Splint (AMS) / Reline	Full Arch AMS (MCI)	Gelb
Material	Methyl Methacrylate	Thermoplastic	Thermoplastic	Methyl Methacrylate + Hydroplastic	Methyl Methacrylate + Hydroplastic	Methyl Methacrylate
Method of manufacture	Heat Cured Process	Thermoform	Thermoform	Thermoform + Acrylic Platform	Thermoform + Acrylic Platform	Acrylic Build Up
Aesthetics	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
Strength	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★	★★★★★
Add-on options <small>(Colours, Clasps, Name and Number Inserts)</small>	✓	✓	✗	✗	✗	✓
Benefits	- Strength (High) - Full Contact - Full Coverage	- Comfort (Soft Lining) - Full Contact - Full Coverage	- Comfort - Full Coverage	- Anterior Only - Comfort - Easy To Adjust Lining	- Full Coverage - Comfort - Easy To Adjust Lining	- Comfort - Day wear - Strength
Recommended use	- Bruxism - Long Term Wear - TMJ Dysfunction	- Bruxism - Introduction To Splint Wearing - Restorative Protection	- Minor Bruxism - Hard Splints are Intolerable (Fit & Comfort)	- Headaches - Bruxism - TMJ Dysfunction	- Headaches - Bruxism - TMJ Dysfunction	- TMJ Dysfunction - Day Wear
Warranty <small>(All 12 Months)</small>	✓	✓	✓	✓	✓	✓

Now Accepting Digital Impressions!



DIGITAL ORTHODONTICS

- The first step to accurate fitting orthodontics is a perfect impression
- Removes steps which multiply inaccuracies



5 lab days turnaround

(09) 887 0327

racedental.co.nz