

## Premium Restorations

Material

Method of manufacture

Aesthetics

Shades / colouring

Strength

Add-on options

(fit to denture, attachments etc. Charges apply)

Benefits

Recommended use

Warranty

			
	Crystalite™	e.max®	Zirconia (Layered)
Material	Monolithic Super Translucent Zirconia	Ivoclar™ e.max® lithium disilicate	Zirconia (layered) coping with porcelain overlay
Method of manufacture	Digitally designed and manufactured	Digitally designed and manufactured	Digitally designed & manufactured (Framework) & built up by hand
Aesthetics	★★★★★	★★★★★	★★★★★
Shades / colouring	No internal effects standard shades	Customised shading	Internal effects customised shading multi-layered
Strength	★★★★★	★★★★★	★★★★★
Add-on options (fit to denture, attachments etc. Charges apply)	Fit to denture	×	Fit to denture, rests and ledges
Benefits	High durability extreme strength metal-free, exceptional accuracy	Premium aesthetics excellent fit metal-free	Good fit metal free
Recommended use	Full mouth, single crowns	Full mouth, crowns, inlays/onlays, veneers, 3/4 crowns	Full mouth, crowns & bridges
Warranty	15 years	5 years	8 years

**We recommend Digital Impressions!**

All intraoral scanners accepted



Fast 5 lab days  
turnaround



100% Australian  
Made

## Standard Restorations

Material

Method of manufacture

Aesthetics

Shades / colouring

Strength

Add-on options

(fit to denture, attachments etc. Charges apply)

Benefits

Recommended use

Warranty

			
	Opalite™	Full Gold Crown	PFM
Material	Monolithic Translucent Zirconia	Non-precious (gold-coloured) or precious (gold) alloy CoCr (silver coloured)	Non-precious or precious alloy coping with porcelain overlay
Method of manufacture	Digitally designed and manufactured	Traditionally-cast	Traditionally-cast & built up
Aesthetics	★★★★★	★☆☆☆☆	★★★★★
Shades / colouring	No internal effects standard shades	Not required	No internal effects standard shades
Strength	★★★★★	★★★★★	★★★★★
Add-on options (fit to denture, attachments etc. Charges apply)	Fit to denture	✓	✓
Benefits	High durability extreme strength metal-free, exceptional accuracy	Exceptional strength	Traditional method choice of alloys
Recommended use	Full mouth, crowns & bridges	Posterior teeth	Full mouth, crowns & bridges
Warranty	15 years	6 years	6 years

**We recommend Digital Impressions!**

All intraoral scanners accepted



Fast 5 lab days  
turnaround



100% Australian  
Made